Delvotest® Fast BL

Fast test for the detection of β-lactam residues in raw cow milk
Insert: January 2018

Introduction

Delvotest® Fast BL is a receptor assay using lateral flow technology that enables the detection of β-Lactam residues of antibiotics in raw cow milk.

Content

- 4 ampoule racks of 25 ampoules each with labelled binding proteins
- 2 containers (labelled with reading indicators) containing 50 test strips each
- 2 sachets of 50 single use pipettes each (250 µl)
- Instructions insert

Caution

- As this test is extremely sensitive to β-lactams any contamination with these substances should be prevented. It is advisable to thoroughly wash and dry hands before starting the test procedure and work in a clean environment.
- Do not use abnormal-looking milk, eg clotted milk, colostrum milk, mastitic milk, etc.
- Do not mix components of the tests from different boxes (ampoules and strips), as it compromises the test's performance.
- Do not use expired product.
- Do not handle the test roughly to avoid damaging the test strips. Rough handling may affect the quality of the test during reading of the results.
- Handle the strip only by touching the green top part (fig 1)
- The temperature of the tested milk should be between 4 -10 °C.
- Do not place the ampoules in the incubator before adding the milk.
- Always keep strips and ampoules under the same storage conditions.

Instructions for use

1. Turn on the Delvotest® incubator (MiniS) and select Program A (Delvotest® BL 1 step) - it will take at least 15 minutes for the incubator to reach the required temperature of 64°C (±2.0°C). Always refer to the Delvotest® Mini S Incubator Manual.
2. Take a strip container and one ampoule rack, out of refrigeration. Allow the container to warm up to room temperature for a period of 10-15 minutes before opening the strip container.
3. Select the required number of ampoules - Detach one or more ampoules, depending on the number of milk samples to be analyzed. Take care that the aluminum foil from the remaining ampoules is not damaged. Clearly label each ampoule with a sample number. Perforate the aluminum foil of the ampoules carefully.
4. Add the milk sample - Milk samples should be well stirred and representative of the milk to be tested. Take a pipette from one of the sachets and squeeze the upper bulb (A) of the pipette tightly. Introduce the pipette tip into the milk sample and release the upper bulb (A). A milk surplus is drawn into the pipette reservoir (B). If there is no milk in the reservoir (B) or air bubbles in the stem (C), squeeze the upper bulb (A) to empty the pipette into the milk sample and repeat pipetting. Squeeze the upper bulb (A) to drop the milk in the ampoule – do not touch the contents of the ampoule with the stem (C) of the pipette. The volume of milk required for the test is 250µl (the amount contained in the pipette stem). For each sample, use a new and clean pipette. Do not work hastily; avoid contaminating the remaining pipettes with milk or your hands. (fig 2)
5. Shake the ampoule until the milk completely dissolves with the reagents. The milk should become pink and contain no residues of the reagent at the bottom of the ampoule. It could take up to 30 seconds to achieve complete dissolution. Not shaking enough might lead to false results.
6. Add the test strip to the ampoule and initiate incubation with the Delvotest® incubator - Take one test strip from the strip container taking into consideration not contaminate the remaining strips and immediately close the container. Hold the test strip from the upper side indicated by the green area, immediately insert one test strip per ampoule insert the ampoule in the incubator and start the incubator timer at 7 minutes incubation at 64°C (±2.0°C). The lower side of the test strip(s) must be inserted in the milk solution. Remove the test strip from the test ampoule after 7 minutes incubation, when the incubator alarm turns on. Gently remove the sample pad (the paper part in contact with the milk) and read the result. (fig 3)
7. Reading - Read immediately within a maximum of 5 minutes after incubation. For visual reading use the strip container label as guidance. If the intensity of the test line (bottom) is less or equal to that of the control line, the sample is positive (contains antibiotics). If the color of the test line is more intense than the control line, the sample is negative. Should the control line not be present or dotted, the test is invalid. Delaying the reading can lead to a false interpretation of the result. Full detection sensitivity and selectivity information is available via your DSM local representatives. The sensitivity of Delvotest® Fast BL for Penicillin G in raw cow milk is 1 to 3 ng/g (ppb). (fig 4)

Additional information

Further information, including sensitivity for specific substances, can be obtained by contacting DSM local representatives. Instrumental reading can be performed with Delvotest® Reader +. Negative and Positive Controls of 4 ppb Pen G (raw cow milk) are available. Contact your local sales representative for more information on acquisition of these items. DSM advises that all samples interpreted as positive should be confirmed (preferably by means of a microbial assay).

Storage

The tests should be stored upright in the original packaging, in the dark at a temperature between 4°C and 8°C and prevented from freezing. The ampoules and strips can remain at room temperature during the testing process, they should however be returned to the advised storage conditions at the end of each working day. The test kits can be shipped at room temperature up to a total time of 4 days.

Shelf Life

The test's shelf life is 12 months from production date. The expiry date is mentioned on the box label.

FOR USA: NOT INTENDED FOR NCIMS TESTING

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